



Yes, you can.®

➤ Invacare® **Action<sup>2</sup><sub>NG</sub>**

### The Next Generation has arrived!

The all-new **Action<sup>2</sup><sub>NG</sub>** has been developed from the tried-and-tested **Action<sup>3</sup>** concept. The specifications meet the basic needs of the customers whilst refusing to compromise on safety or comfort. The next generation **Action<sup>2</sup><sub>NG</sub>** shares the same essential platform as the **Action<sup>3</sup>**, offering a low total weight for easy handling combined with good driving abilities. Common components are cross compatible across the range making the **Action<sup>2</sup><sub>NG</sub>** fit perfectly into the Action family.

NEW



## Light and comfortable

With a weight of just 15.2 kilos, the *Invacare*<sup>®</sup> **Action<sup>2</sup><sub>NG</sub>** fulfills all essential functional requirements. The stable folding frame means the wheelchair is easy to handle, and the lightweight design makes transportation trouble-free. Cross-compatibility with other members of the Action family ensures that the **Action<sup>2</sup><sub>NG</sub>** can be easily adapted to meet individual requirements.



## Quality of life

The tried-and-tested technology of the Action family products allow the customer to go wherever they desire. Cross-compatibility with **Action<sup>3</sup>** and **Action<sup>4</sup>** makes this an ideal, basic, fleet product.

### Features



#### Swing in and out footrest

The 80 degree footrest swings in and out for easy transfers in or out of the chair.



#### Adjustable center of gravity


The newly-designed rear wheel plate can be adjusted by 5 cm.



#### Anti-tippers

Provide extra safety on slopes or uneven ground.

### Technical data




38, 40,5, 43, 45,5, 48 cm




40,5 / 43 cm



48,5, 51 cm



43 cm




20 - 26 cm




35 - 48 cm



80° / 0 - 80°



7° / 0° - 30°




SW+ 20,5 cm



93,5 cm



82 / 84,5 cm



41 cm



16.6 kg self propel  
15.6 kg transit



max. 125 kg



Crash tested  
ISO 7176-19\*

#### Frame colour



Anthracite Black

#### Upholstery



Black nylon

Please notice colors may vary slightly from those displayed above.

\* Our wheelchairs comply above ISO standards and with ISO 10542, and have additional features providing increased levels of occupant security and safety whilst travelling in a motor vehicle. Wheelchair users should however transfer to the vehicle seat and use the vehicle installed restraint system whenever it is feasible.

All information quoted is believed to be correct at time of print. Invacare<sup>®</sup> reserves the right to alter product specifications without prior consultation.



INVACARE<sup>®</sup> Limited  
South Road · Bridgend Industrial Estate · BRIDGEND CF31 3PY  
Tel. +44 (0) 1656 647327 · Fax +44 (0) 1656 649016  
E-mail: uk@invacare.com · [www.invacare.co.uk](http://www.invacare.co.uk)  
Sales Order E-mail: ordersuk@invacare.com

INVACARE<sup>®</sup> Ireland Ltd  
Unit 5 · Seatown Business Campus  
Seatown Road · SWORDS · County Dublin  
Tel. (353) 181 070 84 · Fax (353) 181 070 85  
Email: Ireland@invacare.com · [www.invacare.ie](http://www.invacare.ie)

