

# Rota Stand

*promoting safer patient handling*

Allows a safe, smooth and controlled swivel transfer

Reduces the need to “hold” a patient during a transfer

Encourages an independent and dignified transfer

The Rota Stand has been developed to allow an assisted transfer from one seated position to another. The frame, shin pad, and turn-table combine to give optimum safety to both carer and patient during the transfer.

- **Height adjustable frame to suit all**
- **Foot pedal break to secure turn-table**
- **Central wheels for easy positioning**
- **Angled shin pad for correct standing position**
- **160kg weight limit**



# Rotastand... Solo

The Rotastand... SOLO revolutionizes sit to stand and pivot transfers. The patented mechanism in the base (see picture) allows a sit to stand and pivot transfer to be carried out by a single nurse or carer thus reducing the time it takes to carry out the task as well as halving the cost of having two carers/nurses present during the transfer.

**The nurse/carers does not need to counterbalance the Rotastand...SOLO as the client pulls themselves to a standing position which means the nurse/carers can now assist the client to stand before rotating the client to a new seated position.**



The Central Wheels are a unique innovation because they allow the Rota Stand to be rolled up to the patient or client even after their feet have been positioned on the product.

- Stabilizing bar to allow a one person transfer
- Height adjustable frame to suit all clients
- Central wheels allow easy positioning of the Rotastand... SOLO
- Large Handlebar for a secure grip
- 200Kg weight capacity

