

## The Kelly stand, by Select Healthcare

- Allows nursing and care staff to manoeuvre, toilet and change a person's clothing quickly and without the need to manually lift. so reducing the risk of back injury associated with such tasks.
- To assist people to stand thus promoting rehabilitation and circulation through the body.



The natural motion of standing is achieved by PUSHING DOWN on the arms or a chair surface from which a person is rising.

Traditional standing aids tend to pull the person upwards which can cause painful drag under the person's arms. The unique handlebars allow the person to assist in their own standing by PUSHING DOWN and then giving support in the standing position.

- 160kg weight limit
- Unique handle bars encourage clients to push down in order to stand
- Height and width adjustable shin pad
- Powered raising, lowering and leg opening facility
- Compatible with NON full body weight slings

